

DENZEL
MICROSURGERY

take-off Gewerbepark 65
78579 Neuhausen ob Eck
Germany

Fon +49 (0)7467.911.20
Fax +49 (0)7467.911.22
info@denzel.com
www.denzel.com



NEWSLETTER

Juli-September 2010

DENZEL

MICROSURGERY

Ceramic Line

Denzel Microsurgery offers you our new Ceramic Line it is a new dimension of high quality Instruments.

The Ceramic Line is a special coating which gives your instruments a high improvement in a very simple way. It works in a low-temperature process. It is useful for every standard **stainless steel** and **Titanium** instruments.

Ceramic Line Effects:

- an unique “black is beauty” color
- no light reflections anymore
- Ceramic coating gives an easy cleaning effect because of its very smooth surface
- the surface is 3x harder than standard stainless steel
- anti-allergenic effect
- tissue-friendly und bio-compatible
- high friction resistance
- longer life time
- better corrosion resistance

Let's see some samples....



NEW!!!



Denzel Microsurgery in Titanium

Meet your client's needs and increase your business with “100% German made Titanium Instruments”

Carrying the familiar Denzel quality signature, we are beginning with a specially designed round bodied handle as well as the popular 3 hole handle. All tip patterns are available.

For a custom design please feel free to contact us